

Product datasheet for 4T19-7E7

Cardiac Troponin T (TNNT2) Mouse Monoclonal Antibody [Clone ID: 7E7]

Product data:

Product Type:	Primary Antibodies
Clone Name:	7E7
Applications:	ELISA
Reactivity:	Bovine, Feline, Goat, Human, Mouse, Porcine, Rabbit, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	troponin T2, cardiac type
Database Link:	Entrez Gene 7139 Human
Synonyms:	cTnT, TnTc, TNNT2



[View online »](#)