

Product datasheet for 4T19-7F4

Cardiac Troponin T (TNNT2) Mouse Monoclonal Antibody [Clone ID: 7F4]

Product data:

Product Type:	Primary Antibodies
Clone Name:	7F4
Applications:	ELISA
Reactivity:	Human, Porcine
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	troponin T2, cardiac type
Database Link:	Entrez Gene 7139 Human
Synonyms:	cTnT, TnTc, TNNT2



[View online »](#)