

## Product datasheet for 4T19-9G6

### Cardiac Troponin T (TNNT2) Mouse Monoclonal Antibody [Clone ID: 9G6]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	9G6
Applications:	ELISA, WB
Reactivity:	Mammalian
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	troponin T2, cardiac type
Synonyms:	cTnT, TnTc, TNNT2



[View online »](#)