

Product datasheet for 4T21-267

Cardiac Troponin I (TNNI3) (+ skeletal) Mouse Monoclonal Antibody [Clone ID: 267]

Product data:

Product Type:	Primary Antibodies
Clone Name:	267
Applications:	ELISA, IHC, IP, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	troponin I3, cardiac type
Database Link:	Entrez Gene 7137 Human
Synonyms:	TNNI3, TNNC1



[View online »](#)