

## Product datasheet for **4T21-596**

### Cardiac Troponin I (TNNI3) (+ skeletal) Mouse Monoclonal Antibody [Clone ID: 596]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	596
Applications:	ELISA, IHC, IP, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	troponin I3, cardiac type
Database Link:	<a href="#">Entrez Gene 7137 Human</a>
Synonyms:	TNNI3, TNNC1



[View online »](#)