

Product datasheet for **4T21-801**

Cardiac Troponin I (TNNI3) Mouse Monoclonal Antibody [Clone ID: 801]

Product data:

Product Type:	Primary Antibodies
Clone Name:	801
Applications:	ELISA, IP, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG3
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	troponin I3, cardiac type
Database Link:	Entrez Gene 7137 Human
Synonyms:	TNNI3, TNNC1



[View online »](#)