

## Product datasheet for 4T21-810

### Cardiac Troponin I (TNNI3) Mouse Monoclonal Antibody [Clone ID: 810]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	810
Applications:	ELISA, IP, WB
Reactivity:	Bovine, Canine, Feline, Goat, Human, Porcine, Rabbit
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	troponin I3, cardiac type
Database Link:	<a href="#">Entrez Gene 7137 Human</a>
Synonyms:	TNNI3, TNNC1



[View online »](#)