

Product datasheet for 4T21-916

Cardiac Troponin I (TNNI3) Mouse Monoclonal Antibody [Clone ID: 916]

Product data:

Product Type:	Primary Antibodies
Clone Name:	916
Applications:	ELISA, IP, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG3
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens troponin I type 3 (cardiac) (TNNI3)
Database Link:	Entrez Gene 7137 Human
Synonyms:	TNNI3, TNNC1
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Stem cell - Pluripotency
Protein Pathways:	Cardiac muscle contraction, Dilated cardiomyopathy, Hypertrophic cardiomyopathy (HCM)



[View online »](#)