

Product datasheet for 4T21-C5

Cardiac Troponin I (TNNI3) (+ skeletal) Mouse Monoclonal Antibody [Clone ID: C5]

Product data:

| | |
|----------------|--|
| Product Type: | Primary Antibodies |
| Clone Name: | C5 |
| Applications: | ELISA, IHC, IP, WB |
| Reactivity: | Human |
| Host: | Mouse |
| Isotype: | IgG2b |
| Clonality: | Monoclonal |
| Formulation: | State: Purified |
| Gene Name: | troponin I3, cardiac type |
| Database Link: | Entrez Gene 7137 Human |
| Synonyms: | TNNI3, TNNC1 |



[View online »](#)