

Product datasheet for 4T24-6B11

Albumin (ALB) Mouse Monoclonal Antibody [Clone ID: 6B11]

Product data:

Product Type:	Primary Antibodies
Clone Name:	6B11
Applications:	ELISA, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	albumin
Database Link:	Entrez Gene 213 Human
Synonyms:	ALB, BSA, HSA, Serum Albumin



[View online »](#)