

Product datasheet for **4TP15-TPO34**

Thyroid Peroxidase (TPO) Mouse Monoclonal Antibody [Clone ID: TPO34]

Product data:

Product Type:	Primary Antibodies
Clone Name:	TPO34
Applications:	ELISA, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens thyroid peroxidase (TPO), transcript variant 1
Database Link:	Entrez Gene 7173 Human
Synonyms:	MSA; OTTHUMP00000115536; TDH2A; thyroperoxidase; TPX
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Autoimmune thyroid disease, Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Jak-STAT signaling pathway, Metabolic pathways, Tyrosine metabolism



[View online »](#)