

## Product datasheet for 4Tr26-2B6

### CD71 (TFRC) (soluble form) Mouse Monoclonal Antibody [Clone ID: 2B6]

#### Product data:

|                |   |
|----------------|---|
| Product Type:  | Primary Antibodies                        |
| Clone Name:    | 2B6                                       |
| Applications:  | ELISA, WB                                 |
| Reactivity:    | Human                                     |
| Host:          | Mouse                                     |
| Isotype:       | IgG2a                                     |
| Clonality:     | Monoclonal                                |
| Formulation:   | State: Purified                           |
| Gene Name:     | transferrin receptor                      |
| Database Link: | <a href="#">Entrez Gene 7037 Human</a>    |
| Synonyms:      | TfR1, p90, Transferrin receptor protein 1 |



[View online »](#)