

## Product datasheet for 5G4-4G5

### GAPDH Mouse Monoclonal Antibody [Clone ID: 4G5]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	4G5
Applications:	ELISA, IF, WB
Reactivity:	Bovine, Feline, Goat, Human, Mouse, Porcine, Rat
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens glyceraldehyde-3-phosphate dehydrogenase (GAPDH), transcript variant 1
Database Link:	<a href="#">Entrez Gene 14433 Mouse</a> <a href="#">Entrez Gene 24383 Rat</a> <a href="#">Entrez Gene 2597 Human</a>
Synonyms:	GAPD, CDABP0047
Protein Families:	ES Cell Differentiation/IPS
Protein Pathways:	Alzheimer's disease, Glycolysis / Gluconeogenesis, Metabolic pathways



[View online »](#)