

Product datasheet for 75-310

THAP1 Mouse Monoclonal Antibody [Clone ID: N325B/65]


Product data:

Product Type:	Primary Antibodies
Clone Name:	N325B/65
Applications:	IF, WB
Recommend Dilution:	Immunoblot (IB) Immunocytochemistry (ICC)
Reactivity:	Human
Host:	Mouse
Isotype:	IgG3
Clonality:	Monoclonal
Immunogen:	Fusion protein amino acids 1-213 (full-length) of human THAP1 (also known as THAP domain-containing protein 1, accession number Q9NVV9). Mouse: 92% identity (197/213 amino acids identical) Rat: 92% identity (196/213 amino acids identical) <40% identity with THAP2, THAP3 and other THAP proteins
Formulation:	State: Purified
Gene Name:	THAP domain containing 1
Database Link:	Entrez Gene 55145 Human
Synonyms:	4833431A01Rik; DYT6; FLJ10477; MGC33014
Note:	USERS will cite the UC Davis/NIH NeuroMab Facility in any publication(s) describing the research utilizing the MATERIALS. The suggested acknowledgment statement is as follows: "The monoclonal antibody _ was developed by and/or obtained from the UC Davis/NIH NeuroMab Facility, supported by NIH grant U24NS050606 and maintained by the Department of Neurobiology, Physiology and Behavior, College of Biological Sciences, University of California, Davis, CA 95616." Also, please include the complete clone number (e.g., N52A/42) and the Antibody Registry identification number (e.g., RRID:AB_2120479) to avoid ambiguity. View Research License Agreement



[View online »](#)

Product images:



Transfected cell immunoblot: extracts of COS cells transiently transfected with V5/His-tagged THAP1 or untagged Kv2.1 plasmid and probed with N325B/65 TC supe (left) and mouse anti-V5 (right).