

Streptavidin-AP Kit

Ready-to-use Streptavidin Conjugated Alkaline Phosphatase Solution

Storage: 2-8°C

Catalog No: D29-110
D29-18110 ml
18ml **Intended Use:**

This enzyme conjugate solution can be used to label biotinylated antibodies.

Kit Component:

The enzyme conjugate is provided in ready-to-use form in two sizes.

Catalog No.	Size
D29-110	RTU Streptavidin-AP 110ml
D29-18	RTU Streptavidin-AP 18ml

Recommended Protocol:

- 1) Apply two drop (100ul) or enough volume to cover the tissue section completely
- 2) Incubate at room temperature for 10 minutes.
- 3) Resume standard immunohistochemistry staining procedure or find the procedures on our web site.

Storage:

Store at 2-8°C. Return the reagent back to 2-8°C after each use.

Remarks:

For research use only. Not for diagnostic or therapeutic use.