

Streptavidin-HRP Kit

Ready-to-use Streptavidin Conjugated Horse Radish Peroxidase Solution

Storage: 2-8°C

Catalog No:	D25-500	500 ml	<input checked="" type="checkbox"/>
	D25-110	110 ml	<input type="checkbox"/>
	D25-18	18ml	<input type="checkbox"/>

Intended Use:

This enzyme conjugate solution can be used to label biotinylated antibodies.

Kit Component:

The enzyme conjugate is provided in ready-to-use form in two sizes.

Catalog No.	Size
D25-500	RTU Streptavidin-HRP 500ml
D25-110	RTU Streptavidin-HRP 110ml
D25-18	RTU Streptavidin-HRP 18ml

Recommended Protocol:

- 1) Apply two drop (100ul) or enough volume to cover the tissue section completely
- 2) Incubate at room temperature for 10 minutes.
- 3) Resume standard immunohistochemistry staining procedure or find the procedures on our web site.

Storage:

Store at 2-8°C. Return the reagent back to 2-8°C after each use.

Remarks:

For research use only. Not for diagnostic or therapeutic use.