

Product datasheet for **AB0017-100**

RAB1B


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB 1:500-1:5,000
Reactivity:	Canine, Human, Monkey, Mouse, Rat
Host:	Goat
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Purified recombinant peptide derived from within residues 110 aa to the C-terminus of Rab1b produced in E. coli.
Formulation:	Polyclonal antibody supplied as a 100 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide.
Concentration:	3mg/ml
Purification:	epitope affinity purified
Gene Name:	RAB1B, member RAS oncogene family
Database Link:	NP_112243 Entrez Gene 81876 Human
Background:	RAB1, member RAS oncogene family
Synonyms:	RAB1B
Note:	Detects endogenous levels of Rab1b protein by Western blot in whole cell lysates and transfected cells with GFP-Rab1b cds.
Protein Families:	Transcription Factors



[View online »](#)

Product images:



Anti-Rab1b Ab at 1/1,000 dilution; lysates at 50 ug per lane; rabbit polyclonal to goat IgG (HRP) at 1/10,000 dilution;