


Product datasheet for **AB0126-100**

T7

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB 1:500-1:5,000
Host:	Goat
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Purified recombinant peptide MASMTGGQQMG produced in E. coli.
Formulation:	Polyclonal antibody supplied as a 100 μ l (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide.
Concentration:	3mg/ml
Purification:	epitope affinity purified
Background:	MASMTGGQQMG epitope tag
Note:	The antibody recognizes the T7-tag fused to the amino- or carboxy- termini of targeted proteins in transfected cells.

Product images:



Anti-T7 tag Ab at 2, 500 dilution; rabbit polyclonal to goat IgG (HRP) at 1/10,000 dilution;



[View online »](#)