

Product datasheet for **AB0138-100**

CDH11


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB 1:500-1:2,000
Reactivity:	Zebrafish
Host:	Goat
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Purified recombinant peptide fused with GST and produced in E. coli.
Formulation:	Polyclonal antibody supplied as a 100 µl (4 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide.
Concentration:	4mg/ml
Purification:	epitope affinity purified
Gene Name:	cadherin 11
Database Link:	NP_001788 Entrez Gene 1009 Human
Background:	Cadherin 11, type2, OB-cadherin (osteoblast)
Synonyms:	CAD11; cadherin-11; CDHOB; OB; OB-cadherin; OSF-4
Note:	Detects endogenous levels of total CDH11 protein in the whole 3-day embryo lysates by Western blot.
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

Product images:



Endogenous CDH11 detected at 1:2,000 dilution; lysate of 3 day embryo at 50 ug per lane and rabbit polyclonal to goat IgG (HRP) at 1/10,000 dilution.