

Product datasheet for **AM00042PU-N**

EGFR pTyr1173 (incl. pos. control) Mouse Monoclonal Antibody [Clone ID: 9H2]

Product data:


Product Type:	Primary Antibodies
Clone Name:	9H2
Applications:	ELISA, IF, IHC, IP, WB
Recommended Dilution:	ELISA: Use at 0.05 µg/ml. Western Blot: 0.5 µg/ml for HRPO/ECL detection. <i>Recommended buffer:</i> Casein/Tween 20 based blocking and blot incubation buffer. Use cell lysate delivered with this product as Positive Control (See Protocol below). Immunoprecipitation: Use at 1-10 µg per 10 ⁶ vanadate-treated A431 cells. Immunocytochemistry: Use at 1-10 µg/ml. Immunohistochemistry on Frozen and Paraffin Embedded Material
Reactivity:	Canine, Human, Mouse
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Phosphopeptide conjugated to KLH Epitope: phospho-Tyr1173 N A E pY L R V
Specificity:	Clone 9H2 specifically interacts with the 1170 - N A E pY L R V motif corresponding to the major autophosphorylation site of Human EGFR. Antibody AM00042PU-N does not crossreact with the highly homologous pTyr1248 of activated ErbB2.
Formulation:	PBS containing 0.09% Sodium Azide, PEG and Sucrose State: Purified State: Lyophilized purified Ig fraction
Reconstitution Method:	Restore with 1.0 ml H ₂ O for 15 min at RT
Purification:	Purified from Serum-Free Cell Culture Supernatant by Subsequent Thiophilic Adsorption and Size Exclusion Chromatography
Conjugation:	Unconjugated



[View online »](#)

Storage:	<p>For long-term storage, freeze lyophilizate upon arrival (-20°C). Upon reconstitution, aliquote and freeze in liquid nitrogen, reconstituted antibody can be stored frozen at -80°C up to 1 year. Thaw aliquots at 37°C. Thawed aliquots may be stored at 4°C up to 3 months. Avoid repeated freeze / thaw cycles!</p>
Gene Name:	epidermal growth factor receptor
Database Link:	Entrez Gene 1956 Human P00533
Background:	<p>EGFR/erbB receptors are activated upon binding of EGF and EGF-related growth factors such as TGF alpha, beta-cellulin, Hb-EGF, HRG, or NRG. Binding of these ligands leads to receptor homo- and heterodimerization followed by autophosphorylation and activation of downstream signal transduction pathways (MAPK, PI3K/PKB, and STAT). In addition, EGFR becomes fully activated after phosphorylation of Y869 by src family kinases. Phosphorylation of Y1069 leads to association with cbl and subsequent receptor degradation. Phosphorylation of S1071 by CamKinase II leads to attenuation of kinase activity; phosphorylation of T678 (by PKC) and T693 (by MAPK, p38) interferes with receptor endocytosis/recycling.</p>
Synonyms:	Epidermal growth factor receptor, EGF Receptor, erbB-1, c-ErbB-1
Note:	<p>A Positive Control cell lysate is provided with this antibody (for details see Protocol below).</p> <p>Protocol: Provided Positive Control cell lysate: HepG2 Pervanadate Treated -recommended for Immunoblotting</p> <p>Formulation: Lyophilized cell lysate from HepG2 cells. Serum starved cells were treated for 15 min. with pervanadate. Restore by addition of 200 µl H₂O. After complete solubilization add 200 µl 2x SDS-PAGE sample buffer, mix and incubate at 90°C for 5 min. Store in aliquots at -20°C. Avoid repeated freezing and thawing.</p> <p>Application: 20 µl of positive cell lysate correspond to ca. 80.000 cells. Use 20 µl / lane (mini gel) for HRPO/ECL detection of the target proteins.</p> <p>Please note: The lyophilized cell lysate contains SDS and is not recommended for applications with native proteins such as Immunoprecipitation !</p>

Product images:



Phosphospecificity: Whole cell extracts of control (co), EGF stimulated (EGF) or pervanadate treated (VH) A549 tumor cells were applied to SDS-PAGE (ca 20.000 cells per lane) and transferred to a PVDF membrane. The immunoblot was probed with mab EGFR-9H2 (0.5 ug/ ml) for 1h at RT and developed by ECL (exp. time: 30 sec).