

Product datasheet for **AM00199RP-N**

Phosphotyrosine Mouse Monoclonal Antibody [Clone ID: PY20]

Product data:

Product Type:	Primary Antibodies
Clone Name:	PY20
Applications:	FC
Recommended Dilution:	Flow Cytometry.
Reactivity:	All Species
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Specificity:	Recognizes Phosphotyrosine.
Formulation:	PBS containing 0.09% Sodium Azide and Stabilizing agent Label: PE State: Liquid purified IgG fraction Label: Conjugated to R. Phycoerythrin (RPE)
Concentration:	lot specific
Conjugation:	PE
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light.
Stability:	Shelf life: one year from despatch.



[View online »](#)