

## Product datasheet for **AM00738PU-N**

### Epstein Barr Virus / EBV Early Antigen-R Mouse Monoclonal Antibody [Clone ID: BDI502]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	BDI502
Applications:	ELISA, IHC, IP, WB
Recommended Dilution:	ELISA. Western Blot. Immunohistochemistry on paraffin sections. Immunoprecipitation.
Reactivity:	Epstein Barr Virus
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Infected cell lysate
Specificity:	This antibody recognizes Epstein-Barr virus. Specific for the early antigen (restricted) homologue Bcl-2.
Formulation:	0.01M PBS, pH 7.2 containing 0,09% sodium azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	>90% pure. Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	HHV4, HHV-4



[View online »](#)