

Product datasheet for **AM00757PU-N**

FeLV envelope glycoprotein SU gp70 Mouse Monoclonal Antibody [Clone ID: BDI712]

Product data:

Product Type:	Primary Antibodies
Clone Name:	BDI712
Applications:	ELISA, IF
Recommended Dilution:	ELISA. IFA.
Reactivity:	Feline Leukemia Virus
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	FeLV viral lysate
Specificity:	This antibody recognizes FeLV glycoprotein 70 (gp70).
Formulation:	0.01M PBS, pH 7.2 containing 0,09% sodium azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	>90% pure. Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	P06752
Synonyms:	Feline leukemia virus env, gp70, p15E



[View online »](#)