

Product datasheet for **AM00814PU-N**

Astrovirus Type 1 Mouse Monoclonal Antibody [Clone ID: B192M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B192M
Applications:	ELISA, WB
Host:	Mouse
Isotype:	IgG3
Clonality:	Monoclonal
Immunogen:	Inactivated Astrovirus, serotype 1
Specificity:	This antibody is specific for Astrovirus, serotype 1.
Formulation:	PBS, pH 7.2, containing 0.09% Sodium azide State: Purified State: Liquid Ig fraction
Concentration:	lot specific
Purification:	Protein G chromatography (>95% pure by SDS-PAGE)
Conjugation:	Unconjugated
Storage:	Store the antibody at 2 - 8 °C. Centrifuge vial before opening.
Stability:	Shelf life: one year from despatch.
Background:	A genus of small, circular RNA viruses in the family ASTROVIRIDAE. They cause gastroenteritis and are found in the stools of several vertebrates including humans. Transmission is by the fecal-oral route. There are at least eight human serotypes.
Synonyms:	Human astrovirus-1, HAstV-1
Note:	ELISA. Western blot.



[View online »](#)