

Product datasheet for **AM00854SU-N**

Gram Negative Endotoxin (E. coli O:111 B4 J5) Mouse Monoclonal Antibody [Clone ID: 308]

Product data:

Product Type:	Primary Antibodies
Clone Name:	308
Applications:	ELISA
Recommended Dilution:	May be used in immunoassays to detect and quantitate gram negative endotoxin or cells that express endotoxin.
Host:	Mouse
Isotype:	IgG3
Clonality:	Monoclonal
Immunogen:	E. coli O:111 B4 J5 cells
Specificity:	This antibody reacts with E. aerogenes, S. marcescens, P. mirabilis, A. calcoaceticus and Ps. aeruginosa.
Formulation:	State: Ascites State: Liquid Ascites containing 0,09% sodium azide
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)