

Product datasheet for **AM00915PU-N**

E6 Mouse Monoclonal Antibody [Clone ID: BF7]

Product data:

Product Type:	Primary Antibodies
Clone Name:	BF7
Applications:	ELISA, IHC, R, WB
Recommended Dilution:	Immunohistochemistry. Western BLOT. RIA. ELISA.
Reactivity:	Human Papilloma Virus 18
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	HPV18 E6-Beta galactosidase fusion protein
Specificity:	This antibody reacts to HPV 18 early protein (E6).
Formulation:	10mM Phosphate, pH 7.4 containing 150mM Sodium chloride and 0.06% sodium azide State: Purified State: Liquid purified Ig
Concentration:	lot specific
Purification:	>90% pure (SDS-PAGE). Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one week or (in aliquots) at -20°C for longer. If aliquoted for long term storage, fill volume should be equal to or greater than 50% of the nominal fill volume of the vial used. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	Entrez Gene 1489088 P06463
Synonyms:	Human Papilloma Virus Early Protein 6



[View online »](#)