

## Product datasheet for **AM00917PU-N**

### Human IgE Mouse Monoclonal Antibody [Clone ID: 090-11368]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgE Mouse Monoclonal Antibody [Clone ID: 090-11368]
Clone Name:	090-11368
Applications:	ELISA
Recommended Dilution:	<b>ELISA.</b> <i>Recommended pair for Sandwich Immunoassay</i> Capture: AM00919PU-N Detection: AM00917PU-N
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Formulation:	10mM Phosphate, pH 7.4 containing 150mM Sodium Chloride State: Purified State: Liquid purified Ig fraction (>90% pure by SDS-PAGE) Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)