

Product datasheet for **AM00918PU-N**

Human IgE Mouse Monoclonal Antibody [Clone ID: 090-11348]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgE Mouse Monoclonal Antibody [Clone ID: 090-11348]
Clone Name:	090-11348
Applications:	ELISA
Recommended Dilution:	ELISA. <i>Recommended antibody pair for Sandwich Immunoassay</i> Capture / Detection AM00919PU-N / AM00918PU-N
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Formulation:	10mM Phosphate, pH 7.4 containing 150mM Sodium Chloride State: Purified State: Liquid purified Ig fraction (>90% pure by SDS-PAGE) Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)