

Product datasheet for **AM00920PU-N**

Human IgE Mouse Monoclonal Antibody [Clone ID: ME.114]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgE Mouse Monoclonal Antibody [Clone ID: ME.114]
Clone Name:	ME.114
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Formulation:	PBS, pH 7.2 State: Purified State: Liquid purified Ig
Concentration:	lot specific
Purification:	>95% pure. Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.



[View online »](#)