

Product datasheet for **AM00921HR-N**

Human IgE Mouse Monoclonal Antibody [Clone ID: 1A2]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgE Mouse Monoclonal Antibody [Clone ID: 1A2]
Clone Name:	1A2
Applications:	ELISA
Recommended Dilution:	ELISA: 1/500-1/1,000. Can be used as a Detection antibody in Sandwich ELISA (paired with 4C3).
Reactivity:	Human
Host:	Mouse
Immunogen:	Chimeric Human IgE
Isotype:	IgG1
Formulation:	PBS containing 0.01% Thimerosal as preservative Label: HRP State: Liquid purified Ig fraction Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Note:	Use of Sodium Azide as a preservative will substantially inhibit the enzyme activity of Horseradish Peroxidase.



[View online »](#)