

Product datasheet for **AM00921PU-N**

Human IgE Mouse Monoclonal Antibody [Clone ID: 1A2]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgE Mouse Monoclonal Antibody [Clone ID: 1A2]
Clone Name:	1A2
Applications:	ELISA, WB
Recommended Dilution:	ELISA (1 µg/ml). Western blot (1 µg/ml).
Reactivity:	Human
Host:	Mouse
Immunogen:	Mouse-Human chimeric anti-NP IgE.
Isotype:	IgG1
Formulation:	PBS containing 0.01% Thimerosal as a preservative. State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)