

Product datasheet for **AM00923PU-N**

Human IgG4 (Fc specific) Mouse Monoclonal Antibody [Clone ID: BDI541]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG4 (Fc specific) Mouse Monoclonal Antibody [Clone ID: BDI541]
Clone Name:	BDI541
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot. Immunohistochemistry on frozen sections.
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified human myeloma immunoglobulin. Hybridization of Sp2/0-Ag 14 myeloma cells with spleen cells from Balb/c mice
Isotype:	IgG1
Formulation:	PBS, pH 7.2 State: Purified State: Liquid purified Ig
Concentration:	lot specific
Purification:	Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.



[View online »](#)