

Product datasheet for **AM00926PU-N**

Influenza A (Nucleoprotein) Mouse Monoclonal Antibody [Clone ID: BDI375]

Product data:

Product Type:	Primary Antibodies
Clone Name:	BDI375
Applications:	ELISA, IF
Recommended Dilution:	IFA. ELISA.
Reactivity:	Influenza A Virus
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Infected cell lysate
Specificity:	This antibody recognizes Influenza A virus. Specific for Nucleoprotein.
Formulation:	0.01M PBS, pH 7.2 containing 0,09% sodium azide State: Purified State: Liquid purified Ig
Concentration:	lot specific
Purification:	>90% pure. Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Influenza A Virus, Seasonal Flu



[View online »](#)