

## Product datasheet for **AM00928PU-N**

### Influenza A (Nucleoprotein) Mouse Monoclonal Antibody [Clone ID: B252M]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B252M
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Influenza A Virus
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified influenza virus type A strain H1N1. Hybridization of Sp2/0 myeloma cells with spleen cells from BALB/c mice.
Specificity:	This antibody reacts with Nucleoprotein (NP) of Influenza A. Does not cross react with nucleoprotein of Influenza B virus.
Formulation:	PBS, pH 7.4 containing 0,09% sodium azide State: Purified State: Liquid purified Ig
Purification:	>90% pure (SDS-PAGE). Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.
Synonyms:	Influenza A Virus, Seasonal Flu



[View online »](#)