

## Product datasheet for **AM00947PU-N**

### Influenza B (Nucleoprotein) Mouse Monoclonal Antibody [Clone ID: B275M]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B275M
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Western Blot.
Reactivity:	Influenza B Virus
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified influenza virus type B
Specificity:	This antibody reacts to Nucleoprotein (NP) of influenza virus type B.
Formulation:	PBS, pH 7.4 containing 0,09% sodium azide State: Purified State: Liquid purified Ig
Purification:	>90% pure (SDS-PAGE). Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.
Synonyms:	Influenza B Virus



[View online »](#)