

## Product datasheet for **AM00954PU-N**

### Influenza B (Haemagglutinin HA2) Mouse Monoclonal Antibody [Clone ID: B193M]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B193M
Applications:	ELISA, WB
Recommended Dilution:	<b>ELISA.</b> <b>Western Blot.</b>
Reactivity:	Influenza B Virus
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Purified Influenza B virus
Specificity:	This antibody reacts to Influenza B Virus hemagglutinin 2 (HA2 subunit). The antibody has been tested with the following strains: Strain : Results B/Tokyo/53/99 : Positive. B/Leningrad/86/93 : Positive. B/Victoria/504/00 : Positive. A/Shangdong/9/93 : No Reaction. A/New Caledonia/20/99 : No Reaction.
Formulation:	PBS, pH 7.4 State: Purified State: Liquid purified Ig fraction (>90% pure SDS-PAGE) Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Synonyms:	Influenza B Virus



[View online »](#)