

Product datasheet for **AM00955PU-N**

Influenza B (Haemagglutinin HA2) Mouse Monoclonal Antibody [Clone ID: B198M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B198M
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Western Blot.
Reactivity:	Influenza B Virus
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Purified Influenza B virus
Specificity:	This antibody reacts to Influenza B Virus hemagglutinin 2 (HA2 subunit).
Formulation:	PBS, pH 7.4 containing 0,09% sodium azide State: Purified State: Liquid purified Ig
Concentration:	lot specific
Purification:	>90% pure (SDS-PAGE). Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.
Synonyms:	Influenza B Virus



[View online »](#)