

Product datasheet for **AM00957PU-N**

Influenza B (Nucleoprotein) Mouse Monoclonal Antibody [Clone ID: BDI081]

Product data:

Product Type:	Primary Antibodies
Clone Name:	BDI081
Applications:	IF, LF
Recommended Dilution:	IFA.
Reactivity:	Influenza B Virus
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Hong Kong strain CDC
Specificity:	This antibody recognizes Influenza B virus. Specific for a nuclear protein.
Formulation:	0.01M PBS, pH 7.2 containing 0,09% sodium azide State: Purified State: Liquid purified Ig
Concentration:	lot specific
Purification:	>90% pure. Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Influenza B Virus



[View online »](#)