

Product datasheet for **AM01000PU-N**

Norovirus (capsid protein) Mouse Monoclonal Antibody [Clone ID: B114M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B114M
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Norovirus
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified native Norwalk virus, strain 8FIIa
Specificity:	This antibody recognizes the capsid protein of Norovirus (native and recombinant). Does not cross react with Rotavirus, Adenovirus or Astrovirus.
Formulation:	0.01M PBS, pH 7.2 containing 0,09% sodium azide State: Purified State: Liquid purified Ig
Concentration:	lot specific
Purification:	>90% pure. Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)