

Product datasheet for **AM01107PU-N**

Dehydroepiandrosterone-sulphate Mouse Monoclonal Antibody [Clone ID: BGN/1177/1241]

Product data:

Product Type:	Primary Antibodies
Clone Name:	BGN/1177/1241
Applications:	ELISA
Recommended Dilution:	ELISA: 1/90K - 1/360K.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	HS-BSA conjugate.
Specificity:	This antibody reacts to DHEAS (100%), DHEA (13%), Progesterone (0.01%) and Testosterone (0.02%).
Formulation:	PBS buffer, pH 7.2 containing 0.1% Sodium Azide State: Purified State: Liquid purified Ig
Concentration:	lot specific
Purification:	Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	DHEAS, DHEA 3-Sulfate



[View online »](#)