

## Product datasheet for **AM01139PU-N**

### Haemophilus influenzae Type B Mouse Monoclonal Antibody [Clone ID: 1079/83]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	1079/83
Applications:	ELISA, IF
Recommended Dilution:	<b>ELISA.</b> <b>Immunofluorescence.</b>
Reactivity:	Haemophilus influenzae
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Native, NCTC 7279
Specificity:	This antibody is reactive with <i>Haemophilus influenzae</i> B. Reacts with with isolates fixed by Phenol, Formalin and boiling.
Formulation:	PBS State: Purified State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein G
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)