

Product datasheet for **AM01150PU-N**

Human IgA (heavy chain) Mouse Monoclonal Antibody [Clone ID: Mc24-2E11]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgA (heavy chain) Mouse Monoclonal Antibody [Clone ID: Mc24-2E11]
Clone Name:	Mc24-2E11
Applications:	ELISA, IHC, IP, WB
Recommended Dilution:	ELISA. Western Blot. Immunoprecipitation. Immunohistochemistry on Frozen Sections. Immunohistochemistry on Paraffin Sections (10 µg/ml).
Reactivity:	Human
Host:	Mouse
Immunogen:	Reduced and alkylated Human IgA. Spleen cells from immunised BALB/c mice were fused with cells of the mouse NS1 myeloma cell line.
Isotype:	IgG1
Formulation:	PBS, pH 7.4 State: Purified State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein G
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)