

Product datasheet for **AM01151BT-N**

Monkey IgA Mouse Monoclonal Antibody [Clone ID: IgA5-3B]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Monkey IgA Mouse Monoclonal Antibody [Clone ID: IgA5-3B]
Clone Name:	IgA5-3B
Applications:	ELISA
Recommended Dilution:	ELISA: This antibody may be used in a Direct ELISA or as a Detection antibody in a Sandwich ELISA together with product code SM1508P as a Capture reagent.
Reactivity:	Canine, Feline, Human, Monkey
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS, pH 7.2 containing 0.05% Sodium Azide Label: Biotin State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)