

## Product datasheet for **AM01151PU-N**

### Monkey IgA Mouse Monoclonal Antibody [Clone ID: IgA5-3B]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Monkey IgA Mouse Monoclonal Antibody [Clone ID: IgA5-3B]
Clone Name:	IgA5-3B
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Western Blot: Detects a band of approximately 55kDa under reducing conditions.
Reactivity:	Feline, Monkey
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS, pH 7.4 containing 0.05% Sodium Azide State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)