

## Product datasheet for **AM01152AP-N**

### Human IgA1 (alpha chain) Mouse Monoclonal Antibody [Clone ID: B3506B4]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgA1 (alpha chain) Mouse Monoclonal Antibody [Clone ID: B3506B4]
Clone Name:	B3506B4
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry on frozen sections. ELISA. Western Blot. Immunoblotting.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Formulation:	50mM Tris/ 1mM MgCl <sub>2</sub> containing 0.09% Sodium azide and 50% Glycerol Label: AP State: Liquid purified Ig
Conjugation:	AP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)