

Product datasheet for **AM01152HR-N**

Human IgA1 (alpha chain) Mouse Monoclonal Antibody [Clone ID: B3506B4]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgA1 (alpha chain) Mouse Monoclonal Antibody [Clone ID: B3506B4]
Clone Name:	B3506B4
Applications:	ELISA
Recommended Dilution:	ELISA: 1/500-1/2000.
Reactivity:	Human
Host:	Mouse
Immunogen:	Human IgA1 heavy chains.
Isotype:	IgG1
Formulation:	PBS, pH 7.4 containing 50% Glycerol Label: HRP State: Liquid purified Ig Label: Horseradish Peroxidase
Conjugation:	HRP
Storage:	Store the antibody undiluted at -20°C. Avoid repeated freezing and thawing.



[View online »](#)