

Product datasheet for **AM01153BT-N**

Human IgA2 Mouse Monoclonal Antibody [Clone ID: A9604D2]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgA2 Mouse Monoclonal Antibody [Clone ID: A9604D2]
Clone Name:	A9604D2
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Western Blot. Immunoblotting.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: Biotin State: Liquid purified IgG fraction.
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)