

Product datasheet for **AM01153HR-N**

Human IgA2 Mouse Monoclonal Antibody [Clone ID: A9604D2]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgA2 Mouse Monoclonal Antibody [Clone ID: A9604D2]
Clone Name:	A9604D2
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1/500-1/2000.
Reactivity:	Human
Host:	Mouse
Immunogen:	Human IgA2 heavy chains
Isotype:	IgG1
Formulation:	PBS containing 50% Glycerol Label: HRP State: Liquid purified Ig Label: Horseradish Peroxidase
Conjugation:	HRP
Storage:	Store the antibody undiluted at -20°C. Avoid repeated freezing and thawing.



[View online »](#)