

Product datasheet for **AM01154PU-N**

Human IgD Mouse Monoclonal Antibody [Clone ID: IADB6]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgD Mouse Monoclonal Antibody [Clone ID: IADB6]
Clone Name:	IADB6
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Formulation:	Borate buffered saline without preservatives or stabilizers State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at -20°C only. Avoid repeated freezing and thawing.



[View online »](#)