

Product datasheet for **AM01156PU-N**

Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: SB2H2]

Product data:

| | |
|-----------------------|--|
| Product Type: | Secondary Antibodies |
| Product Name: | Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: SB2H2] |
| Clone Name: | SB2H2 |
| Applications: | ELISA |
| Recommended Dilution: | ELISA. |
| Reactivity: | Human, Rabbit |
| Host: | Mouse |
| Isotype: | IgG1 |
| Formulation: | Tris buffered saline State: Purified State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide |
| Concentration: | lot specific |
| Purification: | Affinity Chromatography on Protein A |
| Conjugation: | Unconjugated |
| Storage: | Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. |



[View online »](#)